

Revision no.: 00

Preparation date: May 08, 2026

Product name	Resilience [®] 1000
INCI name	Human Adipose Derived Stem Cell Conditioned Media/1,2-Hexanediol
전성분명	인체지방조직유래줄기세포배양액/1,2-헥산다이올
Packaging Unit	500mL
Packaging material	Bottle: PET, Cap: HDPE
Storage condition	2~8°C

Test items	Specification	Test method
Appearance	Lightly pink transparent liquid	Visual(in house)
Identification(1)	Absorbance should be observed at 595nm in measuring total protein contents	IR(in house)
Identification(2)	EV(Extracellular vesicles) should be observed by NTA analysis.	NTA analysis(in house)
pH	6.5 ~ 8.5	pH meter(in house)
Sterility test	No microbial growth	KP
Mycoplasma testing	Negative	KP
Adventitious Virus Testing	Negative	KP
Purity test(1): Absence Test for Source Cells and Source Tissues	Negative	KP
Purity test(2): Absence Test for Antibiotics and Serum	Negative	KP
Assay(1): VEGF(ng/mL)	more than 10.0	ELISA(in house)
Assay(2): TGF beta(ng/mL)	more than 0.1	ELISA(in house)
Assay(3): HGF(ng/mL)	more than 1.0	ELISA(in house)

Quality Unit Manager



Factory: #1101~1103, 25, Gasan digital 1-ro, Geumcheon-gu, Seoul, Republic of Korea

R&D center: #1101, 1105A, 25, Gasan digital 1-ro, Geumcheon-gu, Seoul, Republic of Korea

Tel.: 82-2-3402-7495

Fax.: 82-504-279-3586